

## VACPRO Portable Suction with rechargeable battery

Highly versatile, light weight and compact VACPRO portable suction unit delivers powerful suction on demand. It can be used in multiple clinical settings, from the Emergency Room, patient transport, EMS and ambulance to the Recovery Room. The sealed lead-acid battery and hospital-grade, dual-mode charging system keep the unit constantly ready for use.



## Model: VACPRO

### Features:

- ▶ Small- Only 11.5" x 9.8" x 9.8"
- ▶ Light- Just 7.3 lbs
- ▶ Powerful: Max 600 mmHg (+/- 3%)
- ▶ Completely self-contained; protects against contamination of motor and gauge
- ▶ Provides immediate high performance suction with 0-600 mmHg vacuum and 20 lpm flow
- ▶ Running time in battery mode is approximately 50 mins
- ▶ Easy-to-grip handle
- ▶ Collection canister: 800 ml

### Standard Accessories:

- |                                |                      |
|--------------------------------|----------------------|
| ▶ Suction canister x 1pc       | ▶ Suction tube x 1pc |
| ▶ Silicone connect tube x 2pcs | ▶ AC recharger x 1pc |
| ▶ Bacteria filter x 1 pc       | ▶ Carry bag x 1 pc   |



*All specifications are subject to change without further notice.*