



Finger Type Pulse Oximeter

OxiSmarter I



OxiSmarter I is a unique designed Finger Type Pulse Oximeter which can check patient's heart rate and blood-oxygen saturation levels simply and quickly. The OxiSmarter I finger type pulse oximeter (weight 50g only) provides reliable SpO2 and pulse rate measurement and can be worn around the neck. It's convenient for spot check readings. The Large LCD display allows for use in a wide range of light and dark environments. The value of SpO2 and PR can be clearly read in the different angle of LCD display. It's great SpO2 diagnostic tool for physician or nurse to go on an inspection around sickrooms of Hospital.

Features

- Large LCD display with SpO2 waveform
- Display direction is selectable
- Automatic power off after removing finger in 8 secs
- Low batter capacity indication
- Light weight less than 50 gs
- Two AAA batteries for easy replacement
- Brightness can be adjusted in 10 level
- Can be worn around the neck for convenient spot checking

Technical Specification

Patient Range

Adult, Pediatric

SpO2

Range 0% - 99%
 Resolution 1%
 Accuracy 70% - 100%: ± 2%
 0% - 69%: unspecified
 Method Dual wavelength LED

Heart (Pulse) Rate

Range 30 - 240 bpm
 Resolution 1 bpm
 Accuracy ± 2bpm or ± 2%

Display

Type LCD Display
 6 display modes
 Parameters SpO2, PR, Pleth bar, Pleth waveform
 Information Low Battery

Environmental

Operating Temperature 5°C - 40°C
 Storage Temperature -10°C - 40°C
 Operating Humidity 15%RH - 95%RH
 Storage Humidity 10%RH - 95%RH

Classification

Type of protection internally powered equipment
 Safety IEC 60601-1


Mechanical

Size 61 x 32 x 37 mm
 Weight 50g (including 2 batteries)

Battery

Type 2 AAA size (1.5Vdc) Alkaline Batteries or 2 AAA (1.5Vdc) rechargeable batteries
 Operation time: minimum 30 hours for normal Operation

Display Modes



Acare Technology Co., Ltd.

1F, No.28, Lane 3, Guangming Rd., Sansia Town, Taipei County, 237, Taiwan

Tel: +886-2-8674-3370

Fax: +886-2-2672-1363

Email: info@acaretech.com

www.acaretech.com